



The
Most
Intelligent
Skin
Analyzer



1010001101011
1010110101000
10100011011
10101000
30 Million
Clinical Database


WiFi
WiFi Sharing


Cloud Sharing



Shooting, Analyzing, Displaying
3 Functions In 1



1500W HD Lens



24 Times
HD Shadowless LED



www.spectronixmedical.com



UNIQUE ADVANTAGES

- ✓ 180° turnover ultra-sensitive touch screen, makes multi-directional display available.
- ✓ 1500W pixel HD lens, makes images more clear and accurate.
- ✓ 10-inch HD touch screen ensures display more clearly.
- ✓ 24 times the energy-saving HD shadowless LED ensures image without distortion.
*With comprehensive suggestion and customized solutions to be more flexible and changeable.
- ✓ 30 million clinical database as reference, bases on large database analysis, combined with artificial intelligence and image analysis technology, ensure analysis more accurate, with analysis speed only 20 seconds ensures time-saving and high efficient.
- ✓ No need to install software, system runs automatically once switch on, no need for software installation support.
- ✓ Product relates to analysis report scientifically, ensuring auto-recommendation more reasonable.
- ✓ Comprehensive management systems of Advertisement, Client, Product and Analysis.
- ✓ WiFi connection, Cloud sharing, Batch management, Data sharing can be achieved.

SPECIFICATIONS

Power: 20W
 Voltage: 85-246 V/AC
 NG/WG: 6.5 KG/9.0 KG
 Equipment Size: 57*48*43 cm
 Packing: Carton Boxes

WHY TRUESCAN ?

TRUESCAN is the world's most advanced skin detection equipment together with shooting, analyzing, displaying 3 functions in 1. It adopts RGB, UV, PL 3 spectrums, combined with artificial intelligence and image analysis technology. With 10 years market test and over 30 million clinical database analysis, TRUESCAN can in 20 seconds efficiently identify 10 dermis and epidermis skin problems, such as: spots, pores, moisture, textures, wrinkles, pigmentation, and so on. According to the test results, TRUESCAN can auto-recommend products and treatment for customized skin management. It is an essential skin detection equipment for cosmetic company, beauty salon, skin clinic training school, beauty center, beauty shop and spa, the market demand of it can be up to 100%

MARKET COMPARISON

Comparison	Skin Analyzer+Computer	Skin Analyzer +IPAD	TRUESCAN
Spectrum	RGB+UV	RGB+UV+PL	RGB+UV +PL
Working mode	Requires additional computer cable connection(customer own computers)	Need to purchase a specific type of IPAD	Together with Professional tablet, no need buy anymore
Camera pixels	800W	1200W/1500W	1500W
System update	No	No	Auto-Updating
Automatic product recommendation	No	Yes	yes
Product information pre-record	Can not be classified	Can be classified	Can be classified,can be qualitative, can be linked
Analysis speed	1.5 mins	1-2mins	20-30 seconds
Sharing way	Traditional printer	Bluetooth sharing	WIFI sharing
Compatibility	Only WIN7 XP	Only specific IPAD	No need consider compatibility issues
Cost performance	★ ★ ★	★ ★ ★ ★	★ ★ ★ ★ ★

DATA SPEAKS LOUDER THAN WORDS

