



THERM V GOLD

Radio Frequency
Feminine/Skin Treatment Device

Recover Your Confidence with Your "Therm V Gold"

VAGINA or **FACE**

Take your confidence back!



For export only

Radio Frequency Therapy



Radio Frequency Therapy

Take your confidence back!

THERMV**GOLD** Provides two features at the same time!

- * Disposable tips: safe and convenient
- * The system automatically recognizes the assembly of the tips for each type of the treatments on the hand piece, ensuring safer and more efficient treatments.
- * Serves the need of the modern ladies who need discrete solutions with quicker results.

FUNCTION

VAGINA THERAPY


01

Radio Frequency Vaginoplasty Treatment Device

Through the thermal stimuli using the radio frequency, this advanced vagina treatment device enhances the vagina condition and improves the stress due to urinary incontinence with superb efficacy.

- Tightening/whitening effect from the inside to the opening of the vagina.
- Restores the elasticity and contraction of the muscles in the vagina.
- Improves dryness of vagina and atrophic vaginitis.
- Improves blood circulation, and vaginitis / Relieves the pain and improves the sexual feelings feelings during sexual intercourse.
- Generation of collagen in the dermis layer and lifting of labium.
- It does not only treat you but also restores your self-confidence at the same time!

Control UI



>> Select the radio frequency, vaginoplasty treatment video



>> Video clip on vaginoplasty treatment




A. Hand Piece : Use safely and separately for the different treatment menus.
 B. Disposable tip : Maintains the clean hygienic status.

02 FACE THERAPY


Radio Frequency SMILE Lifting Treatment Device

The smile lifting treatment does not leave any damage on the surface of the skin. It shoots Radio Frequency waves to the inside of the mouth (on the mucous layer inside the mouth) to enhance the formation of the collagen, which is highly effective not only in removing the wrinkles but also in enhancing the elasticity of the skin and lifting.


- Improvement of the face wrinkles and skin tightening to restore your confidence and beauty.
- Forms collagen in the deep dermis layer of the skin to maintain your skin healthy and elastic.
- Helps amount of the blood circulation increase, which stimulates the skin fibroblasts, collagen, elastin, and hyaluronic acid generation.
- A new-concept skin treatment that satisfies the recent trend of becoming more beautiful secretly



UI



>> Select the radio frequency, SMILE lifting treatment video



>> SMILE lifting treatment video




Specification

Classification	Specification
Size	240 X 330 X 200mm
Weight	3.5KG
Trolley weight	16KG
Frequency	50/60Hz
RF frequency	460 kHz \pm 5kHz
Rated voltage	AC 100/240V
RF current sending method	Mono-polar, RET
Outputting type	Continuous outputting
Output power	50 Watts
Set temperature range	36~45°C
LCD	Embedded 8 Inch

DISTRIBUTOR

Spectronix
spectrum for medical excellence

 Spectronix Trading LLC, 204-B, Falcon Building, Al Mina Road, Dubai, UAE.

 +971 (0) 4 354 2905

 info@spectronixglobal.com

 www.spectronixmedical.com