

INTRODUCING...

THE 4TH DIMENSION IN NON-SURGICAL TREATMENTS

ULTRALIFT^{4D} *Gold*



“ Now See and Feel the Results instantly”


Face Lift . Double Chin Reduction . Body Shaping . Arms & Butt Lifting

NEW ULTRALIFT 4D GOLD

Ultralift 4D Gold is an ideal treatment for those who do not want to undergo a painful invasive surgery or a thread-lift. This technology uses multiple-depth cartridges to deliver HIGH INTENSITY FOCUSED ULTRASOUND on deeper layers of skin and fat. Each point where ultrasound hits is called a thermal zone and there is instant SKIN-CONTRACTION and cavitation effect (fat liberation) on that point. Compared to conventional HIFU technology which uses 1 line in 1 shot, Ultralift 4D Gold uses 11 lines spaced 10mm apart in a single shot. This enhances the performance of this technology 11 times more than any previous

HIFU MACHINE, with more thermal zones precisely spaced 10mm apart there is more uniform skin contraction giving visible results that the patients can see and feel instantly.


1 SHOT = 11 SHOTS




WHAT IS HIFU ?

ULTRASOUND-THERAPY HIFU FACELIFT: HIFU stands for high intensity focused ultrasound and is an innovative, painless and non-invasive way to reduce wrinkles that lasts up to two years and is skin-friendly.

Additional applications include face and brow lifting, tightening of the connective tissue and cellulite treatment. Highly concentrated ultrasound waves penetrate deep into the skin and stimulate collagen and elastin production. 20% of the result can be seen immediately after the treatment. 3-6 months later, the final result is visible. treatment is completely painless.



SINGLE SHOT OF ULTRALIFT 4D



SINGLE SHOT OF NORMAL HIFU


The highly focused ultrasound activates and accelerates the body's own healing process in the face, neck and décolletage areas and thereby helps sagging skin to tighten itself and become more firm. This method requires no creams, no artificial filler products or toxins. It relies solely on the body's natural regenerative capacity. Since the treatment targets deeper layers, the surface of the skin is not affected. After the treatment, you can simply go about your usual activities.

The Ultrasound Therapy procedure stimulates collagen production via a targeted input of ultrasound energy into the deep tissue layers of the skin. This warms the tissue, which in turn triggers a natural skin response that initiates a regenerative process through the production of neo-collagen. Lasers generally target the uppermost layers of the skin (e.g., to counteract fine lines, wrinkles, and changes in pigmentation). The Ultrasound Therapy HIFU procedure bypasses the skin's surface and focuses entirely on the deep, structural layers of the skin where collagen is deposited. The Ultrasound Therapy HIFU

MAIN APPLICATIONS

- Wrinkle reduction
- Non-surgical face-lifting
- Correction of the facial shape
- Lifting of the cheeks
- Eyebrow lifting
- Lifting of frown and forehead lines
- Lifting of nasolabial folds
- Chin reduction
- Lacrimal sac reduction
- Tightening of the décolletage
- Skin-smoothing treatment post liposuction

BEFORE & AFTER



WHAT'S NEW WITH 4D HIFU?

4D HIFU technology is 11 times more efficient than previous generation technology. Introducing multiple lines in each shot not only enhances the results, it also makes the result more long lasting and visible.




The machine comes with 2 cartridges Focused on depth of 3mm and 4.5 mm but you have an option of buying 5 additional cartridges for the body, makes this system even more USEFUL and Efficient, Excellent results are seen on the abdomen, arms, thighs, hips, love handles.






AUTHORISED DISTRIBUTOR



 +971 (0) 4 354 2905

 info@spectronixglobal.com

 www.spectronixmedical.com

 Spectronix Trading LLC, #205-B, Falcon Building, Al Mina Rd
Near Jumeira 1, Bur Dubai, Dubai, United Arab Emirates.