

Cyclone



**FOR NON SURGICAL
BODY SHAPING & FACELIFTING**

Hexa-Powered Fat Unblocking and Drainage Technology

- ⦿ Non surgical body shaping
- ⦿ Skin tightening – body
- ⦿ Skin Lifting - Face, Arms & Hips
- ⦿ Lymphatic detox & drainage for body & face
- ⦿ Face lifting & wrinkle reduction
- ⦿ Powerful drainage of clogged subcutaneous fat



CYCLONE TREATMENTS

WHY CYCLONE?

Over the years different technologies evolved like vacuum massage, RF, cavitation, etc and each one of them showed some improvement on a particular set of patients and didn't show any improvement on another different set of people.

No single technology can promise results in non-surgical body shaping. Cyclone is a next generation smart body contouring technology which utilizes the experience of past and mixes it with some new super-effective ingredients making it an amazing fat churning technology.

HOW DOES CYCLONE BLEND IN TO THE MOST EFFECTIVE NON SURGICAL BODY SHAPING PROTOCOL EVER?

LET'S DIVE DEEP IN TO THE PROTOCOL ...

1) 360 DEGREE ROTATING TWISTER HANDPIECE

This handpiece combines strong suction with radio frequency with a capability of self-rotating itself along the AXIS as we do the treatment on different areas. It has 3 different heads giving the operator an option to treat any area easily.

The combined effect of STRONG SUCTION coupled with radio frequency and self-rotation, delivers deep heat in the fat layer coupled with strong pinching and twisting effect that easily breaks the stubborn fat and cellulite. This increases the vasodilation of the capillaries to enhance easy drainage of toxins and liberated fat.



2) 40 KHZ ULTRASONIC CAVITATION COUPLED WITH RF

The Strong Ultrasonic Waves are further Propelled with the heat of Radio Frequency which in conjunction pass strong Rippling Effect in the Fat layer blasting the fat out of triglycerides.

3) CUPPING RF BODY

The Strong Suction Effect of RF coupled with Cupping Technique further improves the blood Circulation facilitating Super-Draining of Toxins and Fatty acids.

4) MAGNETIC RF

This technology deliver Very DEEP HEAT to stimulate collagen production and Skin tightening.
Deep heat also further enhances the Metabolism and stimulates body to naturally burn fat at rapid Pace.

5) MAGNETIC RF FACE

The technology is similar to the one used on body , It promotes Collagen developments , promotes development of elastin fibres and henceforth you can see effective skin tightening and wrinkle reduction.

6) CUPPING RF FACE

Unclogs all facials vessels promoting Drainage of toxins and improved blood circulation , especially under eye and double chin area.



CLINICAL DATA




TECHNICAL SPECIFICATIONS

Voltage	110v/220v 50HZ/60 Hz
Power	800 W
Rotation	360°
Vacuum	0-90 Kpa
Ultrasonic	40 KHz
RF	5 MHz
Depth	0.2 MM - 3 MM
Handle	6 ↑
Package size	50.5x52x113 cm
Gross weight	60 Kg



AUTHORISED DISTRIBUTOR



 +971 (0) 4 354 2905

 info@spectronixglobal.com

 www.spectronixmedical.com

 Spectronix Trading LLC, #205-B, Falcon Building, Al Mina Rd
Near Jumeira 1, Bur Dubai, Dubai, United Arab Emirates.